SoClean. & SoClean 2 Go Compatibility Chart

Most CPAP machines with a standard hose will work with the SoClean right out of the box. Some CPAP machines may need a specific adapter, particularly units with a heated hose. To see whether you need to have a special adapter to use the SoClean with your CPAP, please refer to our guide below.



N/A = Not Applicable
✓ = Out of Box
+ = Need Adapter
X = Not Compatible

		PAP	WITH HUMIDIFIER	WITH HEATED HOSE		
Aeiomed						
0 00	Everest	•	~	N/A		
Apex Medical						
500	ICH	•	~	N/A		
	XT	•	•	N/A		
Carefusion						
	PureSom	•	~	N/A		
Devilbiss						
	ComfortPAP	•	~	N/A		
Fisher & Paykel						
	Icon	•	~			
	SleepStyle 200 Series	~	~	N/A		
	SleepStyle 600 Series	+	+			
Philips Respironics						
	PR System One REMstar 60 Series	•	~			



	N/A = Not Applicable	✓ = Out of Box	x + = Need Adapte	er X = Not Compatible
		PAP	WITH HUMIDIFIER	WITH HEATED HOSE
	PR System One REMstar			
		~	✓	N/A
	SleepEasy	✓	✓	N/A
	SleepEasy II	✓	✓	N/A
	M-Series V 2	~	~	N/A
	M-Series V 3	✓	✓	N/A
	Remstar	✓	✓	N/A
PMI Probasic				
	Zzz-PAP	•	✓	N/A
Puritan Bennet				
	GoodKnight 420	•	✓	N/A
	Sandman	•	✓	N/A
RESmart				
	600	✓	✓	N/A
	700	•	✓	N/A
ResMed				
	S8	•	✓	N/A
Par Je	S9	~	~	
	AirSense 10	+	+	
Transcend				
	Transcend II	✓	(HME) X	N/A

